00173 Rome – Via Giacomo Delitala, 106 Italy
Tel. (+39)06/88640509 – Fax (+39)06/88384287
VAT: 09673041001
Website www.aktive-corp.com

INSTRUCTIONS FOR USE (IFU)

SUMMARY

BACKFLUSH.....	3
CHANDELIER.....	5
CHARLES CANNULAE.....	7
FIBER OPTIC, ILLUMINATED PIC.....	9
HIGH PRESSURE INJ/ASP SYSTEM.....	11
INFUSION ASPIRATION TUBING SET.....	13
INFUSION CANNULA AND ANTERIOR MANTAINER.....	15
IRIS RETRACTORS.....	17
KNIFES_DELAMINATION SPATULA_RETINAL SCRAPER.....	19
LASER FIBER ILLUMINATED.....	20
LASER FIBER_LASER FIBER EXTENSIBLE.....	22
NEEDLES AND CANNULAE.....	24
PHACO NEEDLE.....	26
RETINAL TIP.....	28
TROCAR.....	30
TWEEZER.....	32
VITRECTOMY CUTTER.....	34

Aktive®

Ophthalmic Instruments

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Rev. 5 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION ACTIVE/PASSIVE BACKFLUSH

EN

With fixed or extendable terminal

INDICATIONS

The “ACTIVE/PASSIVE BACKFLUSH” is a sterile single-use medical device intended for ophthalmic surgery. The ACTIVE/PASSIVE BACKFLUSH SYSTEM models are available in variable gauges of 20g, 23g, 25g, 27g, 29g, in the following types

ACTIVE/PASSIVE BACKFLUSH SYSTEM Standard terminal	ACTIVE/PASSIVE BACKFLUSH SYSTEM Protected terminal	ACTIVE/PASSIVE BACKFLUSH SYSTEM Brush terminal
	Fixed terminal Extendable terminal	Fixed terminal Extendable terminal

BACKFLUSH devices in all types and models are intended exclusively for the transfer of fluids from the inside to the outside of the eye and occasionally in reverse during Vitreo-Retinal surgery. Therefore, fluids (either aspiration or infusion) will pass through the BACKFLUSH MD during Vitreo-Retinal surgery. To be used by specialised personnel only.

SAFETY WARNINGS

Please note that the BACKFLUSH MD must only be used by ophthalmic surgeons, who are the only ones who are familiar with their characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.

- A) Use the device before the indicated expiry date, only for its intended use.
B) Do not use the product if the packaging is not intact or inadvertently opened before use.
C) Open aseptically.
D) Store in a dry place away from heat, dusty or humid areas. The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or re-use.
Any reprocessing procedures (decontamination washing and sterilisation)

cause irreparable degradation of the MD materials exposing the operator and patient to serious safety risks. In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and to the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. The pack contains 1 active/passive BACKFLUSH handpiece complete with M/F connector and cap. Depending on requirements, the BACKFLUSH can be connected directly to the mechanical suction of the Vitrectomy equipment if the machine connector is Female. Otherwise, use the Male/Female converter to convert the machine connector and make the connection to the BACKFLUSH possible. The BACKFLUSH handpiece is terminated with a Male Luer Cone connector.

This makes it possible to use the BACKFLUSH in an ACTIVE configuration, transferring the suction/infusion parameters of the vitrectomy machine directly to the device terminal. By occluding/opening the silicon support through the finger, the surgeon will be able to divert the aspiration/infusion alternately from the terminal of the Backflush to the hole of the silicon support, effectively stopping or starting the action at the tip of the backflush.

Use in PASSIVE configuration is possible by occluding the handpiece with the supplied LUER CAP. In this method of use, the BACKFLUSH will only act as a Passive outflow route to liquids inserted by other routes. Occlusion or disocclusion of the silicone holder by the surgeon will allow the infusion to flow or be blocked. Using the system passively opens up the possibility that during the period of occlusion of the silicone holder the intraocular pressure may rise due to the pressure of the fluid that other systems are infusing. It is therefore important for the surgeon to assess this aspect when the infusion (external to the MD and not controlled by it) is left active and the Backflush MD is used with the silicone support occluded.

Depending on the model, the backflush terminal could be standard, or fitted with a cylindrical silicone tube- or a brush-type terminal. Both the cylindrical and brush-type silicone terminal can be of the fixed or extendable type. In the latter case, use the slider on the handpiece to extend the terminal to the desired length. It is important for the surgeon to assess the density of the fluid to be aspirated or infused in relation to the gauge used. This is based on the operator's experience. In case of incorrect assessment, if the fluid is too dense, occlusion of the device itself may occur.

DISPOSAL

Used devices should be disposed of according to standard hospital procedures to control the risk of contamination and taking care to protect operators from unintentional injury during disposal. Use rigid collection buckets for used and contaminated sharp MDs. Secondary packaging, IFUs, blisters and Tyvek® pouches, are medical waste that can be assimilated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with municipal sorting rules.

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Rev. 5 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION CHANDELIER LIGHTING

EN

INDICATIONS

The "Chandelier Lighting" (CHANDELIER) is a sterile single-use medical device for ophthalmic surgery. THE FIBRES according to the different models are made in 20g, 23g, 25g, 27g gauges with single or split terminals. Types of CHANDELIER:

single CHANDELIER	double CHANDELIER
The optical fibre with CHANDELIER terminal is single	The optical fibre with CHANDELIER terminal is double
For each model, connectors are available in each gauge according to different connection standards	
TYPE: A-B-T-D-S-R-K-G-E	
The chandelier can be supplied in combination with a TROCAR	

Regardless of the model or type, the device is indicated for the transfer of illumination from an external light source (vitrectomy instrument) to which the fibre connects, to the posterior pole of the eyeball during vitreoretinal surgery.

SAFETY WARNINGS

Please note that the CHANDELIER MD should only be used by ophthalmic surgeons, who are the only ones who are familiar with its characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements performed by the doctor/surgeon during the procedure.

- A. Use the device before the indicated expiry date, only for its intended use.
B. Do not use the product if the packaging is not intact or inadvertently opened before use.
C. Open aseptically.
D. Store in a dry place away from heat, dusty or humid areas.
E. The CHANDELIER is a passive means of transferring the light signal emitted by external machinery. Exposure times and illumination powers are to be referred to the manufacturer of the instrument.

The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or re-use.

Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing both operator and patient to serious safety risks.

In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. The package contains 1 CHANDELIER LIGHTING.

Connect the CHANDELIER LIGHTING to the supplied instrument using the prepared connector. Ensure that the CHANDELIER LIGHTING is fitted with the appropriate connector for the machine in use. Check the lighting output on the terminal. The terminal has a rubber ring near it, which acts as a stop during insertion.

If the surgeon wishes to have a dedicated insertion route, a pre-assembled TROCAR can be included in the package.

The TROCAR is already pre-inserted on a handpiece equipped with a cutting blade so that the surgeon can perform the Trocar incision/insertion manoeuvre in a single movement. The insertion technique of the Trocar is variable according to the different surgical techniques, but we recommend an insertion movement tangent to the tissue in the initial part of the insertion and then verticalising the angle after a few mm. This type of insertion if performed correctly (please refer to various surgical techniques to be consulted by the practitioner) will make the use of suture at the end of the operation unnecessary. As soon as the incision has been made with the handpiece with blade and TROCAR pre-inserted, perform the reverse movement of the handpiece+blade. This movement will result in the TROCAR being positioned in situ through the incision. Take care to carry out the movement slowly but progressively, possibly holding the TROCAR now inserted in the Sclera with a fine-bladed forceps.

The TROCAR is equipped with a valve, which can be removed if necessary by the surgeon, should its removal be necessary during the operation. Avoid any contact of the blade and Trocar with elements of the surgical area, even if sterile. The extremities are very delicate and could hardly withstand even the slightest impact without deforming and losing their effectiveness.

Insert the fibre inside the eyeball through the previously made incision or through the previously positioned Trocar in the case of minimally invasive surgery. Gently penetrate into the bulb, taking care to move the retaining ring to the correct distance from the retinal surface in order to provide diffuse light without hindering surgical manoeuvres. It is important for the operator to carefully check the portion of the CHANDELIER LIGHTING inserted, acting on the retainer, to avoid making contact with the retina.

In case of DOUBLE CHANDELIER LIGHTING, the procedure to be performed is duplicated, although the connector remains single.

It is recommended that the CHANDELIER LIGHTING MD in the chosen calibre be combined with Trocar of the same calibre, which is complementary to it. In fact, the use of the MD with a smaller gauge (larger gauge) Trocar will result in entry obstruction. On the contrary, the use of Trocar with a larger gauge (smaller gauge), although it has no contraindications for the MD, turns out to be technically incorrect for the purposes of surgery, as it will create undesirable outflow pathways with poor pressure holding.

DISPOSAL

The used device must be disposed of according to standard hospital procedures to control the risk of contamination and with the utmost care to protect operators from unintentional injury during disposal. Secondary packaging, IFUs, blister packs and Tyvek® pouches, are medical waste that can be treated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with urban sorting rules.

PACKAGING AND STERILISATION

Single, sterile pack, packaged in aseptically openable blister packs.

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Rev. 5 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION CHARLES FLUTE NEEDLES

EN

INDICATIONS

The ‘CHARLES FLUTE NEEDLES’ are a sterile single-use medical device intended for ophthalmic surgery. CHARLES FLUTE NEEDLES models are available in gauges varying from 20g, 23g, 25g, 27g, in the following types

CHARLES FLUTE NEEDLES Standard terminal	CHARLES FLUTE NEEDLES Protected terminal	CHARLES FLUTE NEEDLES Brush terminal
---	--	--

“CHARLES FLUTE NEEDLES” MDs in all types and models are exclusively intended for the transfer of fluids from the inside to the outside of the eye and occasionally in reverse during Vitreo-Retinal surgery.

Therefore, fluid (both in aspiration and infusion) passes through the “CHARLES FLUTE NEEDLE” MD during Vitreo-Retinal surgery. To be used by specialised personnel only.

SAFETY WARNINGS

Please note that the “CHARLES FLUTE NEEDLE” MD must only be used by ophthalmic surgeons, who are the only ones who are familiar with its characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.

- A) Use the device before the indicated expiry date, only for its intended use.
B) Do not use the product if the packaging is not intact or inadvertently opened before use.
C) Open aseptically.
D) Store in a dry place away from heat, dusty or humid areas. The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or reuse. Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing the operator and patient to serious safety risks. It is important for the surgeon to assess the density of the fluid to be aspirated or infused in relation to the gauge used. This is based on the operator’s experience. In the event of incorrect

assessment, if the fluid is too dense, occlusion of the device itself may occur.

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METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. The package contains 1 CHARLES FLUTE NEEDLE with Luer connector. The CHARLES FLUTE NEEDLE is a device to be used in conjunction with the BACKFLUSH handpiece, of which it is the terminal. Connect the Luer carefully and make sure the connection is perfect and complete.

Depending on the model, the CHARLES FLUTE NEEDLE may be standard or fitted with a cylindrical silicone or brush-type terminal tube.

DISPOSAL

The used device must be disposed of in accordance with standard hospital procedures to control the risk of contamination and taking care to protect operators from inadvertent injury during disposal. Use rigid buckets to collect used and contaminated sharp MDs. The secondary packaging, IFUs, blisters and Tyvek® pouches, are medical waste that can be assimilated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with municipal sorting rules.

PACKAGING AND STERILISATION

Single, sterile pack, packaged in aseptically openable blister packs. ETO Sterilisation

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Rev. 5 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION OPTICAL FIBRE FOR LIGHTING standard – wide angle – shielded illuminated pic

EN

INDICATIONS

The “OPTICAL FIBRE FOR LIGHTING” is a sterile, single-use medical device for ophthalmic surgery. The FIBRE according to the different models are available in gauges varying from 20g, 23g, 25g, 27G, with variable tip terminals. Types of OPTICAL FIBRES:

Standard Fibre The fibre tip is truncated and polished flush with the metal end in which it is contained	Wide Angle Fibre The fibre tip protrudes from the end containing it in a semi-elliptical shape	Shielded Fibre The semi-elliptical shaped fibre tip partially protrudes from the tube containing it for only a portion, the tube being cut at 45°.
For each model, connectors are available in each gauge according to different connection standards TYPE : A-B-T-D-S-R-K-G-E		
The fibre may be fitted with a spherical end cap to act as an illuminated scleral depressor.		
The optical fibre may be fitted with a 45° cutter terminal on the guide tube terminal. The MD will be called ILLUMINATED PIC		

Regardless of the model or type, the indication of the device is the transfer of illumination from an external light source (vitrectomy instrument) to which the fibre is connected, to the posterior pole of the eyeball during vitreoretinal surgery. In the case of the ILLUMINATED PIC, the transfer of light is combined with the possibility of being able to incise the retinal tissue.

SAFETY WARNINGS

Please note that FIBER MDs must only be used by ophthalmic surgeons, who are the only ones who know their characteristics and fields of action; any other use may compromise the safety of both user and patient. In addition, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.

- a) Use the device before the indicated expiry date, only for its intended use.
- b) Do not use the product if the packaging is not intact or inadvertently opened before use.
- c) Open aseptically.
- d) Store in a dry place away from heat, dusty or humid areas.

e) Optical fibre is a passive medium for transferring the light signal emitted by external machinery. Exposure times and illumination powers are to be referred to the manufacturer of the instrument.

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METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. The package contains 1 OPTICAL FIBRE FOR LIGHTING or its ILLUMINATED PIC variant, equipped with plastic protection. The protection, which covers the tip of the FIBRE, must be removed before use. Connect the FIBRE to the instrument using the connector provided. Ensure that the fibre is fitted with the appropriate connector for the machine in use. Check the lighting output on the fibre terminal. Insert the fibre into the eyeball through the previously made incision or through the previously positioned Trocar in case of minimally invasive surgery. Penetrate into the bulb gently, taking care not to reach the retinal surface. The fibre is of sufficient length to compensate even particularly myopic eyes. For this reason it is important that the operator carefully checks the inserted portion of the fibre to avoid making contact with the retina.

Depending on the model chosen, the optic fibre can be fitted with a removable spherical end cap, which acts as an illuminated scleral depressor in the event that the operator wants to check the exact correspondence outside the bulb of a retinal break.

The optic fibre may be fitted with a 45° angled cutting terminal at the end of the guide tube (ILLUMINATED PIC denomination). The operator will therefore have a cutting terminal with associated light. In this case, the device will allow the operator to incise the retinal surface and at the same time have precise illumination of the incision point.

It is recommended that the OPTICAL FIBRE FOR LIGHTING MD in the chosen calibre be associated with Trocar of the same calibre, complementary to it. In fact, the use of MD with a Trocar of a lower calibre (larger gauge) will result in entry obstruction. On the contrary, the use of a Trocar with a larger calibre (smaller gauge), although it has no contraindications for the MD, is technically incorrect for the purposes of surgery, as it will create undesirable outflow pathways with poor pressure holding.

DISPOSAL

The used device must be disposed of according to standard hospital procedures to control the risk of contamination and taking care to protect operators from unintentional injury during disposal. For PICs, use rigid collection buckets for used and contaminated sharp MDs. Secondary packaging, IFUs, blister packs and Tyvek® pouches are medical waste that can be assimilated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with municipal sorting rules.

PACKAGING AND STERILISATION

Single, sterile pack, packaged in aseptically openable blister packs. ETO sterilisation

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Rev. 1 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION EN

High pressure infusion/aspiration system

INDICATIONS

The “High Pressure Infusion/Aspiration System” is a sterile single-use medical device for ophthalmic surgery.
The “High Pressure Infusion/Expiration System” models are intended exclusively to provide an automated means of performing 2 types of functions with the aid of external equipment:
- High-pressure infusion of pre-filled viscous fluid into a syringe.
- Suction of viscous fluid from the eye cavity to the outside.
The High Pressure Infusion/Aspiration System is equipped with different types of connectors, which can be adapted to all common air connections of commercially available equipment. The different connectors have the designation A, B, C, D, E, F, G, H.
The device is indicated for infusion and/or aspiration of silicone oil in vitreoretinal surgery. It is supplied in two sizes: for syringes with a plunger diameter of 14.0 mm (AK-0383/X2) and 14.5 mm (AK-0383/X). The 14.5 mm version (AK-0383/X) is also available on request complete with syringe to be filled (AK-0383/XS). The X stands for the selected connector, which, depending on the model, will be represented by the letter A, B, C, D, E, F, G, H. The high-pressure infusion/aspiration system is designed and tested to withstand operating pressure not exceeding 4 BAR. To be used by specialised personnel only.

SAFETY WARNINGS

Please note that the High Pressure Infusion/Expiration System should only be used by ophthalmic surgeons, who are the only ones who are familiar with their characteristics and fields of action; any other use may compromise the safety of the user and patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.
A. Use the device before the indicated expiry date, only for its intended use.
B. Do not use the product if the packaging is not intact or inadvertently opened before use.
C. Open aseptically.
D. Store in a dry place away from heat, dusty or humid areas.
E. Check that the high-pressure infusion/aspiration system connector is suitable for the machine in use.
F. Check that the diameter of the clamping cylinder has been chosen to match the syringe in use (if not using AKTIVE syringe)
G. Do not exceed the operating pressure of 4 BAR
The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or re-use.
Any reprocessing procedures (decontaminating washing and sterilisation)

cause irreparable degradation of the MD materials exposing both operator and patient to serious safety risks.
In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. The pack contains 1 high-pressure infusion/aspiration system.

Method of use in the infusion configuration:

1. Ensure that the pre-loaded syringe is of a suitable diameter to fit the pneumatic cylinder, insert the coupling system with an interlocking movement, bringing the O-ring inside the syringe cylinder. Rotate the coupling system 90° in order to tighten the syringe to the system. Connect the air tube to the air outlet of the machine.
2. Before use, carefully check that you have followed the previous steps and that all connections are stable.

The syringe is ready to infuse automatically, with external machine control. It will be up to the surgeon to use different types of infusion routes (not provided by the system) connected to the syringe via a traditional Luer Lock.

Method of use in the aspiration configuration:

The system must be used with the syringe supplied with AKTIVE.

1. Prepare the empty syringe, without cap, by inserting the plunger (piston + rubber head) to the end of the stroke. Before doing so, wet all components with (sterile) saline solution to facilitate sliding.
2. Unscrew the piston, leaving the rubber head inserted at the bottom of the syringe. Proceed to assemble the High Pressure Infusion/Aspiration System according to step 1. Connect the system tube to the aspiration path of the external machine. The system is now ready for aspiration work. The system can be used to aspirate liquids from the eye cavity. It will be up to the surgeon to use different types of extraction routes (not provided by the system) connected to the syringe via a conventional Luer Lock. In the event that the liquid to be extracted exceeds the volume contained in the syringe, the operator will have to disassemble the syringe and empty it, acting on the rubber head after reassembling the piston. The operator can now proceed to reassemble the system having a second aspiration stroke equal to the volume of the syringe. All operating parameters are dependent on machinery external to the system.
3. Before proceeding with use, carefully check that you have followed the previous steps and that all connections are stable.

DISPOSAL

The used device must be disposed of according to standard hospital procedures to control the risk of contamination and taking care to protect operators from unintentional injury during disposal.
Secondary packaging, IFUs, blister packs and Tyvek® pouches, are medical waste that can be assimilated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with municipal sorting rules.

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Rev. 1 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION INFUSION/ASPIRATION LINE COMBINATION OF SEVERAL LINES IN AN INFUSION/ASPIRATION SYSTEM

EN

INDICATIONS

The “INFUSION LINES” are sterile single-use medical devices for ophthalmic surgery. The models of INFUSION/ASPIRATION LINES consist of tubes in variable gauges and lengths, designed to transport (infuse) or evacuate (aspirate) fluids, saline solution, air, organic residues through connection to bags/flasks or external instruments or both. The combination of several types of INFUSION/ASPIRATION LINES can form an INFUSION/ASPIRATION SYSTEM. The various models are equipped with connectors which, depending on the type and combination, can be of the male or female Luer type. INFUSION/ASPIRATION LINES can be associated with ‘millipore’ type filters for increased liquid or air filtration, directional flow taps, deflectors and cannulae. They are intended exclusively to establish a connection between the outside and inside of the eye in order to ensure the passage of fluids with the support of a machine that controls the dynamics (external device). They can be used for both anterior and posterior segment surgery. To be used by specialised personnel only.

SAFETY WARNINGS

Please note that INFUSION/ASPIRATION LINES should only be used by ophthalmic surgeons, who are the only ones who know their characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/ surgeon performs during the procedure.

- A) Use the device before the indicated expiry date, only for its intended use.
- B) Do not use the product if the packaging is not intact or inadvertently opened before use.
- C) Open aseptically.
- D) Store in a dry place away from heat, dusty or humid areas.
- E) Be certain of the intended connection, the function of the line, the correctness of the connectors, the direction of infusion/ aspiration as well as all flow, flow rate, pressure parameters that will be transferred to the line.
- F) INFUSION/ASPIRATION lines or their combination in the INFUSION/ASPIRATION SYSTEM are a non-specific medical device for simple fluid/air transfer.
- G) The maximum pressure must not exceed 20 psi (1.3 bar), the design limit of the device, for any reason.
- H) Use for pressure-driven, controlled or manual infusion/aspiration and connection to surgical instrumentation is the responsibility of the user.

The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or re-use. Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing both operator and patient to serious safety risks. In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. The package contains one or more INFUSION/ASPIRATION lines or their combination in an INFUSION/ASPIRATION SYSTEM. Connect the Luer-type connectors carefully and make sure the connection is perfect and complete. As the connection is non-specific, it is left to the user, under their own responsibility, to use the INFUSION/ASPIRATION lines or their combination in the INFUSION/ASPIRATION SYSTEM according to the combination considered most useful for the surgical practice used. Remember never to exceed a pressure of 20 psi (1.3 bar). It is important for the surgeon to assess the density of the fluid to be aspirated or infused. This is based on the operator's experience. In case of incorrect assessment, occlusion of the device itself may occur if the fluid is too dense.

DISPOSAL

The used device must be disposed of according to standard hospital procedures to control the risk of contamination and taking care to protect operators from unintentional injury during disposal. Secondary packaging, IFUs, blister packs and Tyvek® pouches are medical waste that can be treated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with urban sorting rules.

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Rev. 5 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION INFUSION LINES EN (Low Pressure Infusion Lines and Chamber Maintainers)

INDICATIONS

The “INFUSION LINES” are sterile single-use medical devices for ophthalmic surgery.

The models of INFUSION LINES are available in gauges of 20g, 23g, 25g, 27g.

They are intended exclusively to establish a connection between the outside and inside of the eye in order to ensure fluid passage.

The construction principle is to ensure the passage of fluids through a flexible external connection pathway (connection pipe). Depending on the fluid and its density, the infusion line can be of the standard or reinforced type.

The types are:

- Infusion line for low-pressure solution (non-reinforced): suitable for the passage of low-density fluids;
- Reinforced infusion line for low pressure solution: designed to guarantee the passage of high density fluids;
- Infusion line for Chamber Maintainer: designed to guarantee the passage of low-density fluids (saline solution);

Can be used for both anterior and posterior segment surgery.

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SAFETY WARNINGS

We would like to remind you that INFUSION LINES must only be used by ophthalmic surgeons, who are the only ones who are familiar with their characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to verify that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.

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- E) Low-pressure (unreinforced) cannulae should only be used for infusion/aspiration of low-density fluids, particularly saline solution or fluids of the same density. Do not use for infusion of denser fluids. The risk is line breakage or disconnection of the tube from the cannula.

Instead, it is recommended that only reinforced cannulae be used if infusion of Silicone Oil or any other fluid denser than water is required.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving table or mother table aseptically. The package contains an Infusion Line or Chamber Maintainer. Attach the Luer carefully and ensure perfect and complete connection. They can be used to infuse saline or silicone oil into the anterior chamber, posterior chamber or both, respectively. They can be inserted, depending on the surgical technique either through the sclera, through a trocar (previously prepared) or directly into the clear cornea through a previously prepared incision.

The saline solution INFUSION LINE and anterior chamber maintainer are intended exclusively for the infusion of lowpressure, low-density fluids.

It is important that the surgeon assesses the density of the fluid to be aspirated or infused in relation to the gauge used. This is based on the operator's experience. In the event of incorrect assessment, if the fluid is too dense, the device itself could become occluded and in the worst cases the line could rupture or the tube could become disconnected from the cannula.

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Secondary packaging, IFUs, blister packs and Tyvek® closures are medical waste that can be treated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with urban sorting rules.

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Rev. 5 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION IRIS RETRACTORS

EN

INDICATIONS

The “IRIS RETRACTORS “(RETRACTORS)” are a sterile single-use medical device for ophthalmic surgery. IRIS RETRACTORS, are intended exclusively for use as a mechanical means of pupil enlargement, in the case of hypertrophy of the iris sphincter or poor action of the iris dilator muscle during ophthalmic surgery. To be used by specialised personnel only.

SAFETY WARNINGS

Please note that IRIS RETRACTORS should only be used by ophthalmic surgeons, who are the only ones who are familiar with their characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements performed by the doctor/surgeon during the procedure. A) Use the device before the indicated expiry date, only for its intended use. B) Do not use the product if the packaging is not intact or inadvertently opened before use. C) Open aseptically. D) Store in a dry place away from heat, dusty or humid areas. The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or re-use. Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing the operator and the patient to serious safety risks. In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, checking the integrity of the package. The package contains 4 (or 5, depending on the packaging chosen) pre-assembled PMMA IRIS hooks with sliding silicone ring. Place the contents on the serving or parent table aseptically. The package may contain 4 to 5 IRIS RETRACTORS. The surgeon must make a number of incisions equal to the number of retractors to be used. The incisions should be made along the circumference of the cornea closest to the limbus. The incision may be of minimum diameter, 18-23G, made with a STAB, Side Port or MVR type cutter. Insert the hooked side of the device inside the wound until it reaches the pupillary sphincter. Perform a ‘mechanical hook’ manoeuvre of the retractor in order to be able to use the retractor as a mechanical means of traction. Take care to leave the pre-mounted silicone ring on the outside of the bulb on the opposite side of the hook and use it as a “stopper” for the desired dilation adjustment. It is customary to use 4 or 5 hooks equally distributed over the 360° corneal circumference. Once the traction of all retractors has been adjusted, the pupil will be mechanically dilated and it will be easier for the surgeon to access the posterior chamber of the eyeball.

DISPOSAL

The used device must be disposed of according to standard hospital procedures to control the risk of contamination and taking care to protect operators from inadvertent injury during disposal. Use rigid collection buckets for used and contaminated sharp MDs. Secondary packaging, IFUs, blister packs and Tyvek® pouches are medical waste that can be treated as municipal waste and therefore disposed of according to procedures applied by the hospital and according to urban sorting rules.

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Rev. 5 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION OPTICAL FIBRE FOR LASERS With illumination

EN

INDICATIONS

The “OPTICAL FIBRE FOR LASER WITH ILLUMINATION” (ILLUMINATED LASER FIBRE) is a sterile single-use medical device for ophthalmic surgery.

The ILLUMINATED LASER FIBRE according to the different models are available in gauges of 20g, 23g, 25g Types

Connectors to the laser source are available in each gauge according to different connection standards TYPE: SMA - QTL -LMN -HGM- DRC

For each model, connectors to the light source are available in each gauge according to different connection standards TYPE: A-B-T-D-S-R-K-G-E

Regardless of the model, the device is indicated for the transfer of laser light and illumination light simultaneously within the posterior pole of the eyeball during vitreoretinal surgery, for the performance of laser photocoagulation, with the aid of illumination that improves the view of the point of impact. Both lights originate from a dual external source (both laser and illumination) to which the ILLUMINATED LASER FIBRE is connected.

SAFETY WARNINGS

Please note that the ILLUMINATED LASER FIBRE should only be used by ophthalmic surgeons, who are the only ones who are familiar with their characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.

A) Use the device before the indicated expiry date, only for its intended use.

B) Do not use the product if the packaging is not intact or inadvertently opened before use.

C) Open aseptically.

D) Store in a dry place away from heat, dusty or humid areas.

E) ILLUMINATED LASER FIBRE is a passive means of transferring the dual laser+light signal emitted by external machinery. Exposure times and powers are to be referred to the instrument manufacturer.

The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or re-use. Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing the operator and patient to serious safety risks.

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METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. The package contains 1 ILLUMINATED LASER FIBRE with plastic protection. The protection, which covers the tip of the ILLUMINATED LASER FIBRE, must be removed before use. Connect the ILLUMINATED LASER FIBRE to the supplied instrument using both prepared connectors (laser and illumination connector). Ensure that the ILLUMINATED LASER FIBRE is fitted with the appropriate connectors for the machine in use. Check the output of the “pilot” illumination (aiming beam) on the terminal of the ILLUMINATED LASER FIBRE by lowering the brightness of the light source. Take care to adjust the light source to a power that still allows the pilot light of the laser (aiming beam) to be displayed, in order to always show the point of impact of the laser if triggered. Insert the ILLUMINATED LASER FIBRE inside the eyeball through the incision previously made or through the Trocar previously positioned in the case of minimally invasive surgery. Penetrate into the bulb gently, taking care not to reach the retinal carpet. The ILLUMINATED LASER FIBRE is of sufficient length to compensate even particularly myopic eyes. For this reason it is important that the operator carefully checks the distance of the terminal to the retinal surface in order to avoid contact.

It is recommended that the OPTICAL FIBRE FOR LASER WITH ILLUMINATION MD in the chosen calibre be combined with Trocar of the same calibre, complementary to it. In fact, the use of the MD with a Trocar of a lower calibre (larger gauge) will result in entry obstruction. On the other hand, the use of a Trocar with a larger gauge (smaller gauge), although it has no contraindications for the MD, is technically incorrect for the purposes of surgery, as it will create undesirable outflow paths with poor pressure holding.

DISPOSAL

The used device must be disposed of according to standard hospital procedures to control the risk of contamination and with the utmost care to protect operators from unintentional injury during disposal. Secondary packaging, IFUs, blister and Tyvek® pouches, are medical waste that can be treated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with urban sorting rules.

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Rev. 5 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION STANDARD AND EXTENSIBLE

EN

Optical fibre for lasers

INDICATIONS

The “OPTICAL FIBRE FOR LASERS” (LASER FIBRE) is a sterile single use medical device for ophthalmic surgery. The LASER FIBRE according to the different models are available in gauges varying from 20g, 23g, 25g, 27G. Types of LASER FIBRES

Standard LASER FIBRE The tip of the LASER FIBRE is polished flush with the metal end in which it is contained	Extensible LASER FIBRE The tip of the LASER FIBRE protrudes from the end containing it with an operator- adjustable curvature
For each model, connectors according to different connection standards are available in each gauge TYPE: SMA - QTL -LMN -HGM- DRC	

Regardless of the model or type, the indication of the device is the transfer of laser light from an external source to which the LASER FIBRE connects to the posterior pole of the eyeball during vitreoretinal surgery, for the performance of laser photocoagulation.

SAFETY WARNINGS

Please note that the LASER FIBRE MD should only be used by ophthalmic surgeons, who are the only ones who are familiar with its characteristics and fields of action; any other use may compromise the safety of both user and patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.

- A. Use the device before the indicated expiry date, only for its intended use.
B. Do not use the product if the packaging is not intact or inadvertently opened before use.
C. Open aseptically.
D. Store in a dry place away from heat, dusty or humid areas.
E. The LASER FIBRE is a passive means of transferring the light signal emitted by external machinery. Exposure times and powers are to be referred to the manufacturer of the instrument.

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Any reprocessing procedures (decontamination washing and

sterilisation) cause irreparable degradation of the MD materials exposing both operator and patient to serious safety risks. In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. The package contains 1 LASER FIBRE with plastic protection. The protection, which covers the tip of the LASER FIBRE, must be removed before use. Connect the LASER FIBRE to the supplied instrument using the supplied connector. Ensure that the LASER FIBRE is fitted with the appropriate connector for the machine in use. Check the output of the “pilot” illumination (aiming beam) on the terminal of the LASER FIBRE. Insert the LASER FIBRE inside the eyeball through the incision previously made or through the Trocar previously positioned in the case of minimally invasive surgery. Penetrate the bulb gently, taking care not to reach the retinal surface. The LASER FIBRE is of sufficient length to compensate even particularly myopic eyes. For this reason it is important that the operator carefully checks the distance of the terminal with the retinal surface in order to avoid contact.

In the case of the extendable LASER FIBRE, the handpiece will be equipped with a cursor. By acting on it, the operator will be able to allow the LASER FIBRE to escape from the metal terminal which acts as a guide. The LASER FIBRE, previously bent with a “memory effect”, will perform a bending movement as it exits, which will cause it to assume a variable angle with respect to the axis of the guide tube. The greater the portion of fibre that exits, the greater the angle of curvature. This condition will allow the operator to orient the LASER beam at an angle to the axis of orientation of the guide tube, performing photocoagulation treatments in hard-to-reach areas.

It is recommended that the Optical Fibre for Laser MD in the chosen calibre be combined with Trocar of the same calibre, complementary to it. In fact, the use of the MD with a Trocar of a lower calibre (larger gauge) will result in entry obstruction. On the contrary, the use of a Trocar with a larger gauge (smaller gauge), although it has no contraindications for the MD, turns out to be technically incorrect for the purposes of surgery, as it will create undesirable outflow pathways with poor pressure holding.

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Rev. 1 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION NEEDLES AND CANNULAE

EN

INDICATIONS

Sterile single-use medical devices included in the large category of "NEEDLES AND CANNULAE" for ophthalmic surgery, constitute access routes to transfer liquids, substances or solids dissolved in fluids from inside the eye to the outside or *vice versa*. They can be a simple vehicle of liquid into a previously created access route (in this case we speak of a cannula) or they can create one before starting the transfer (in this case we speak of a needle).

The only difference is their ability respectively not to be or to be sharp. For all of them, the principle is always the same: to convey liquids, substances or solids dissolved in fluids to or from the eye. Gauge diameters can vary from 19 G up to 38 G. The types and functions are listed below:

	Indication	Area of use
Subretinal injection needle	- Infusion of drugs into the sub-retinal area - Aspiration of sub-retinal oedemas	Vitreoretinal surgery
Bimanual cannulae	Pair of cannulae: - One, for aspiration, intended for aspiration of post-extraction crystalline lens residue - Another, for infusion of aqueous humour replacement fluid, in complementary action to that of aspiration Both cannulae constitute, when connected to external machinery, an irrigation/aspiration system	Cataract surgery
Dual Bore Cannula	Cannula indicated for fluid exchange. While the main infusion route injects fluid into the eyeball (usually perfluorine), the peripheral route delivers the fluid to be replaced (usually saline)	Vitreoretinal surgery
Viscous fluid cannula	Cannula suitable for injecting or aspirating very high density fluids (silicone oil) into the eyeball with the help of a syringe or machine-connected infusion/aspiration system.	Vitreoretinal surgery
Needles and Cannulae: - For anaesthesia - For capsulotomy - For hydrodissection - For irrigation - For Viscoclastic materials - For Air Infusion	To deliver anaesthesia drugs to intra- or peri-ocular areas To incise the anterior capsule of the crystalline lens and infuse saline solution To infuse liquid into the intracapsular space of the crystalline lens To infuse liquid into the anterior chamber To infuse liquid into the anterior chamber To infuse air into the anterior chamber	Cataract surgery
Tear Cannulae	To infuse liquid into the tear duct	Tear duct irrigation

To be used by specialised personnel only.

SAFETY WARNINGS

Please note that the "NEEDLES and CANNULAE" MDs must only be used by ophthalmic surgeons, who are the only ones who are familiar with their characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.

- A) Use the device before the indicated expiry date, only for its intended use.
B) Do not use the product if the packaging is not intact or inadvertently opened before use.
C) Open aseptically.
D) Store in a dry place away from heat, dusty or humid areas.

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Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing the operator and patient to serious safety risks.

In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and to the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. Connect the Luer carefully and make sure the connection is perfect and complete. The use of individual devices responds to precise surgical techniques that cannot be listed in these instructions and remains a personal legacy of the operator's experience and manual dexterity.

It is important that the surgeon assesses the density of the fluid to be aspirated or infused in relation to the gauge used. This is a result of the operator's experience. In case of incorrect assessment, if the fluid is too dense, occlusion of the device itself may occur.

DISPOSAL

Used devices should be disposed of according to standard hospital procedures to control the risk of contamination and taking care to protect operators from unintentional injury during disposal. Use rigid collection buckets for used and contaminated sharp MDs.

The secondary packaging, IFUs, blisters and Tyvek® pouches, are medical waste that can be assimilated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with municipal sorting rules.

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Rev. 7 vom 20.10.2022

MEDICAL DEVICE DESCRIPTION PHACOEMULSIFICATION TIP

EN

INDICATIONS

The “PHACOEMULSIFICATION TIP” (Phaco Tip) is a sterile single-use medical device for ophthalmic surgery. The Phaco Tips according to the different models are made in gauges of 19g, 20g, 21g, 22g, variable tip angles of 0°, 15°, 30°, 45°, length between 22mm and 30mm, metric or withworth threads, star or square spanner. Regardless of the type, the device is indicated for the crushing of the Crystalline through the elongation of the material by a train of ultrasonic (U/S) pulses emitted by a Phaco handpiece to which the tip is connected. Phaco tip types:

Straight-Type Phaco Tip	Curved-Type Phaco Tip (Kellmann)	Curved-Type Phaco Tip
-------------------------	----------------------------------	-----------------------

PHACO TIP devices, in all types and models, are intended exclusively for the mechanical crushing of the natural crystalline material of the eye and its aspiration during cataract surgery. The PHACO TIP MD will be passed through by ultrasounds, which will cause the elongation of the material at the frequency established by the manufacturer of the handpiece to which it is connected; this elongation process is associated with the passage of liquids and surgical residues (in both aspiration and infusion) through the use of the SLEEVE (supplied) during cataract surgery. To be used by specialised personnel only.

SAFETY WARNINGS

Please note that PHACO TIP MDs should only be used by ophthalmic surgeons, who are the only ones who are familiar with their characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.

- Ensure correct choice of model and thread compatibility: do not force.
 - Use the device before the indicated expiry date, only for its intended use.
 - Do not use the product if the packaging is not intact or inadvertently opened before use.
 - Open aseptically.
 - Store in a dry place away from heat, dusty or humid areas.
 - The device is for single use only. Do not re-use or re-sterilise
- The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or re-use. Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing the operator and the patient to serious safety risks. In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. The package contains 1 PHACO TIP with plastic protection, a TEST CHAMBER and a SLEEVE. The protection, which covers the PHACO TIP also serves as a tightening key. Without removing the tip from the guard, thread the device onto the PHACO handpiece, taking care not to force the operation in the initial connection phase. Be sure to use the correct thread (metric or WithWorth) corresponding to the handpiece being used. Gently thread the tip and only tighten it firmly when it has reached its natural stop. It is only at this point that the guard should be separated from the PHACO TIP, which is now screwed onto the handpiece, taking care not to damage the tip during this operation. Use the supplied sleeve to cover the PHACO TIP. This in turn will slide over the tip and must also be screwed onto the PHACO handpiece. Take care not to engage the sleeve's own side holes during assembly. The end of the PHACO TIP should come out of the sleeve with the side holes free and undamaged. It will be the operator's responsibility to fine-tune the last screwing step of the sleeve to adjust the position of the holes and end of the sleeve relative to the tip. Use the TEST CHAMBER to coat the handpiece+tip+sleeve creating a “closed chamber” so as to allow the functional verification and “tuning” of the tip to the handpiece according to the handpiece manufacturer's instructions. Always use the sleeve: non-use would result in lack of irrigation flow, the incision edges would be subjected to harmful mechanical and thermal stress.

DISPOSAL

The used device must be disposed of according to standard hospital procedures to control the risk of contamination and taking care to protect operators from inadvertent injury during disposal. The secondary packaging, IFUs, blisters and Tyvek® pouches, are medical waste that can be assimilated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with municipal sorting rules.

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EN

INDICATIONS

The “RETINAL TIP” is a sterile single-use medical device for ophthalmic surgery. RETINAL TIP models are divided into grasping devices and cutting devices, in gauges varying from 20g, 23g, 25g according to the following table:

RETINAL TIP with gripper at the tip	RETINAL TIP with serrated gripper	RETINAL TIP with straight scissor	RETINAL TIP with curved scissor
Gripping device with gripping surface at the top of the gripper	Gripping device with serrated gripping surface	Straight-type cutting device, crossed blades	Curved-type cutting device, crossed blades

The RETINAL TIP, in all types and models, is intended exclusively for the gripping and/or cutting of retinal tissue and underlying layers during Vitreo-Retinal surgery. The Device is for single use only and must be associated exclusively with the use of a specific multi-use handle that can operate it (AKTIVE code AK-1642). Use on any tissue other than retinal tissue or the layers bordering it, due to its greater consistency, could deform the device or block it in its action. Do not use with handles other than the one indicated. To be used by specialised personnel only.

SAFETY WARNINGS

Please note that the RETINAL TIP should only be used by ophthalmic surgeons, who are the only ones who are familiar with its characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as the instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure. A) Use the device before the indicated expiry date, only for its intended use. B) Do not use the product if the packaging is not intact or inadvertently opened before use. C) Open aseptically. D) Store in a dry place away from heat, dusty or humid areas. The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or re-use. Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing

the operator and the patient to serious safety risks. In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. Take care to protect the tip from any contact with elements of the surgical area, even if sterile. The gripping/cutting tips are very delicate and could hardly withstand even the slightest impact without becoming deformed and losing efficiency. Assemble the RETINAL TIP aseptically to the AKTIVE AK-1642 Multipurpose Handle, which will have been previously sterilised for use on patients. The RETINAL TIP will be attached to the handpiece using the threaded joint ring provided with the handle. Depending on the surgical technique the RETINAL TIP to reach the retinal area of use can be used through an incision made by the surgeon on the eyeball. This incision can be made with a scalpel blade and therefore the RETINAL TIP will pass through the scleral tissue to reach the posterior pole in direct contact with the sclera itself. Differently, in the case of minimally invasive surgery using a Trocar, the RETINAL TIP will be introduced directly through the trans-scleral trocar. It is critical that the RETINAL TIP, whether forceps or scissors, be inserted through the wound or trocar in a “closed” manner so as not to encounter any stopping or friction surfaces, providing the smallest possible diameter for passage. It is recommended that the RETINAL TIP MD in the chosen calibre be combined with Trocars of the same calibre, complementary to it. In fact, the use of RETINAL TIP MD with Trocar of a smaller gauge (larger gauge) will result in the obstruction of entry. On the contrary, the use of Trocar with a larger gauge (smaller gauge), even though it has no contraindications for the MD, turns out to be technically incorrect for the purposes of surgery, as it will create undesirable outflow pathways with poor pressure seal.

DISPOSAL

Used devices should be disposed of according to standard hospital procedures to control the risk of contamination and taking care to protect operators from unintentional injury during disposal. Use rigid collection buckets for used and contaminated sharp MDs. Secondary packaging, IFUs, blister packs and Tyvek® pouches are medical waste that can be treated as municipal waste and therefore disposed of according to procedures applied by the hospital and according to urban sorting rules.

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Rev. 5 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION TROKAR

EN

INDICATIONS

The TROCAR is a sterile single-use medical device for ophthalmic surgery. TROCAR models are available in variable gauges of 20g, 23g, 25g, 27g. TROCAR devices, in all types, are exclusively intended to establish a connection pathway between the outside and the inside of the eye in order to ensure a protected passage of surgical instruments, needles, cannulae and the transfer of fluids from the outside to the inside of the eye or *vice versa* during Vitreo-Retinal surgery. To be used by specialised personnel only.

SAFETY WARNINGS

Please note that TROCAR MDs should only be used by ophthalmic surgeons, who are the only ones who are familiar with their characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.
A) Use the device before the indicated expiry date, only for its intended use.
B) Do not use the product if the packaging is not intact or inadvertently opened before use.
C) Open aseptically.
D) Store in a dry place away from heat, dusty or humid areas.
The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or re-use. Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing the operator and the patient to serious safety risks.
In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical

table aseptically. The pack can contain 1 to 3 TROCARs plus an infusion line. The TROCAR is already pre-inserted on a handpiece equipped with a cutting blade so that the surgeon can perform the Trocar incision/insertion manoeuvre in a single movement. The insertion technique of the Trocar is variable according to the different surgical schools, but we recommend an insertion movement tangent to the tissue in the initial part of the insertion and then verticalising the angle after a few mm. This type of insertion if performed correctly (please refer to various surgical techniques to be consulted by the practitioner) will make the use of suture at the end of the operation unnecessary. As soon as the incision has been made with the handpiece with blade and TROCAR pre-mounted, perform the reverse movement of the handpiece+blade. This movement will result in the TROCAR being positioned in situ through the incision. Take care to carry out the movement slowly but progressively, possibly holding the TROCAR now inserted in the Sclera with a fine-bladed forceps.
The TROCAR is equipped with a valve, which can be removed if necessary by the surgeon, should its removal be necessary during the operation. Avoid any contact of the blade and Trocar with elements of the operating area, even if sterile. The extremities are very delicate and could hardly withstand even the slightest impact without deforming and losing their effectiveness. An infusion line, to be used in connection with the surgical machine being used, is also supplied in the TROCAR package. The infusion line can be inserted into any of the available TROCARs, depending on the most convenient positioning for the surgical phase.
The infusion line has a conical cross-section and therefore it is sufficient to slightly force it into the Trocar (cylindrical crosssection) in order to have sufficient clamping for the passage of saline solution without leakage. The infusion line is intended exclusively for saline solution and is not designed for pressures greater than 1000 mmHG. Use a different infusion line with adequate reinforcement if the Trocar is used to pass fluids denser than saline solution or for pressures greater than 1000 mmHG.
It is important for the surgeon to assess the density of the fluid to be aspirated or infused in relation to the gauge used. This is based on the operator's experience. In the event of incorrect assessment, if the fluid is too dense, occlusion of the device itself and, in the worst cases, rupture of the line or disconnection of the tube from the cannulae may occur.

DISPOSAL

Used devices must be disposed of according to standard hospital procedures to control the risk of contamination and taking care to protect operators from unintentional injury during disposal. Use rigid collection buckets for used and contaminated sharp MDs.
Secondary packaging, IFUs, blister packs and Tyvek® pouches are medical waste that can be treated as municipal waste and therefore disposed of according to procedures applied by the hospital and according to urban sorting rules.

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Rev. 5 vom 10.11.2022

MEDICAL DEVICE DESCRIPTION TWEEZER

EN

INDICATIONS

The “TWEEZER” is a sterile single-use medical device for ophthalmic surgery.

TWEEZER models are divided into gripping devices and cutting devices, in gauges varying from 20g,23g,25g,27g according to the following table:

FOREIGN BODY TWEEZER	GRIP TWEEZER	JAWS TWEEZER
Gripping device for foreign bodies, with 3branch gripper	Gripping device with gripping surface at the top of the gripper	Gripping device with serrated gripping surface
straight CUTTING TWEEZER	curved CUTTING TWEEZER	VERTICAL TWEEZER
Straight cutting device with crossed blades	Curved cutting device with crossed blades	Cutting device with parallel sliding blades of the 'guillotine' type

TWEEZER devices, in all types and models, are intended exclusively for gripping and/or cutting retinal tissue and underlying layers during Vitreo-Retinal surgery. The only exception to this specific purpose is the FOREIGN BODY TWEEZER, designed with 3 branches (instead of 2) to be able to better clamp and extract foreign bodies that have penetrated the eyeball due to various accidents. To be used by specialised personnel only.

Use on any tissue other than retinal tissue or the layers bordering it, given its greater consistency, could deform the device or block it in its action.

SAFETY WARNINGS

Please note that the TWEEZER DMs should only be used by ophthalmic surgeons, who are the only ones who are familiar with their characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as the instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.

- A) Use the device before the indicated expiry date, only for its intended use.

- B) Do not use the product if the packaging is not intact or inadvertently opened before use.

- C) Open aseptically.

- D) Store in a dry place away from heat, dusty or humid areas.

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Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing the operator and the patient to serious safety risks.

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METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or parent table aseptically. Take care to protect the tip from any contact with elements of the surgical area, even if sterile. The gripping/cutting tips are very delicate and could hardly withstand even the slightest impact without becoming deformed and losing efficiency.

Depending on the surgical technique, the TWEEZER MD to reach the retinal area of use can be used through an incision made by the surgeon on the eyeball. This incision can be made with a scalpel blade and therefore the TWEEZER MD will pass through the scleral tissue to reach the posterior pole in direct contact with the sclera itself. Differently, in the case of minimally invasive surgery using a Trocar, TWEEZER will be introduced directly through the trans-scleral trocar. It is essential that TWEEZER MDs, whether Tweezers or scissors, are inserted through the wound or trocar in a “closed” manner so as not to encounter any stopping or friction surfaces, offering the smallest possible diameter for passage.

It is recommended that the TWEEZER MD in the chosen calibre be combined with Trocars of the same calibre that are complementary to it. In fact, the use of the TWEEZER MD with Trocar of a smaller gauge (larger gauge) will result in entry obstruction. On the other hand, the use of Trocar with a larger gauge (smaller gauge), although it has no contraindications for the MD, is technically incorrect for the purposes of surgery, as it will create undesirable outflow pathways with poor pressure holding.

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Secondary packaging, IFUs, blister packs and Tyvek® pouches are medical waste that can be treated as municipal waste and therefore disposed of according to procedures applied by the hospital and according to urban sorting rules.

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Rev. 9 vom 20.10.2022

MEDICAL DEVICE DESCRIPTION VITRECTOMY CUTTER

EN

INDICATIONS

The "VITRECTOMY CUTTER" is a sterile single-use medical device intended for ophthalmic surgery.

The VITRECTOMY CUTTER models are available in gauges varying from 20g, 23g, 25g, 27g. They are divided into 2 types:

- VITRECTOMY CUTTER with active pneumatic action and spring return (1 pneumatic and 1 aspiration tube).
- VITRECTOMY CUTTER with both pneumatic cutting action and return (3 tubes: 2 pneumatics and 1 aspiration).

Both can be equipped with a blade with 1 cutting profile (single cut) or 2 cutting profiles (dual cut). The mechanical cutting action is the same for both types except for a different construction solution, which in the first case employs 1 pneumatic drive path associated with a passive return path (spring) while in the second it employs 2 active pneumatic paths alternately opposed to each other to perform the oscillating blade movement function. In the first case only one tube is needed to operate the blade, in the second 2. The VITRECTOMY CUTTER, regardless of type and model, is intended exclusively for the aspiration of the Vitreous Humor, through a cutting and aspiration action during vitreoretinal surgery, for its partial or total removal. Since the Vitreous, due to the collagen fibres that constitute it, is an elastic substance with high viscosity, it needs to be cut into small portions to be removed. The VITRECTOMY CUTTER performs this function through the use of a blade inside a port at the end of the handpiece with a pneumatically operated guillotine movement associated with aspiration. The guillotine blade placed at the end of the VITRECTOMY CUTTER can be supplied with 1 cutting profile or with double profile, depending on the surgeon's choice. In this second case we will have the designation VITRECTOMY DUAL CUT. The VITRECTOMY CUTTER must be connected to external machinery, from which it is controlled. The number of strokes and aspiration power are determined by the machine based on the surgeon's settings.

WARNINGS

Please note that the VITRECTOMY CUTTER must only be used by ophthalmic surgeons, who are the only ones who are familiar with their characteristics and fields of action; any other use may compromise the safety of the user and the patient. Furthermore, as these instruments are for manual use, it is recommended to check that in the operating area (operating theatre), there are no disturbing elements that could interfere with the movements that the doctor/surgeon performs during the procedure.

- Store in a dry place away from heat, dusty or humid areas.
- Use the device before the indicated expiry date, only for its intended use.
- Do not use the product if the packaging is not intact or inadvertently opened before use.
- Open aseptically.
- Respect the connection logic and do not force the connectors
- Cutting speed and aspiration are to be referred to the manufacturer of the instrument to which it is connected, check compatibility on the basis of the enclosed table.
- Before performing surgery, connect the device and ensure its correct functioning also by visually checking the smooth movement of the guillotine blade.
- Use on any substance or tissue other than Vitreous Humor may impair the functionality of the device.
- If the tip bends or malfunctions, do not proceed and replace the device.
- During use be careful around the retina: it may be subjected to unwanted stresses
- The device is for single use only, do not re-use or re-sterilise
- In order for the VITRECTOMY CUTTER to be perfectly efficient, it is important that the user, before use, ensures to have chosen a model of VITRECTOMY CUTTER in the AKTIVE production line that meets the following characteristics:
 - Connectors of a type complementary to the connections provided by the manufacturer of the instrument to which it is to be connected.
 - That the VITRECTOMY CUTTER, has the CPM (cut per minute) and operating pressure (PSI) parameters equal to the values that the manufacturer declares as the machine's operating standard.
- If both points of letter L are not met, the VITRECTOMY CUTTER will not work properly.
- Before use, carefully check that all previous points have been met

The device is for single use only. The manufacturer does not guarantee the product in case of re-sterilisation or re-use. Any reprocessing procedures (decontamination washing and sterilisation) cause irreparable degradation of the MD materials exposing the operator and patient to serious safety risks. In the event of a serious accident occurring due to the use of the device, the event must be promptly reported to Aktive and to the health authorities.

METHOD OF USE

Open the package and remove the contents in a sterile manner, verifying the integrity of the package. Place the contents on the serving or surgical table aseptically. Take care to protect the tip from any contact with elements of the operating area, even if sterile. The tips are very delicate and could hardly withstand even a slight impact without deforming and losing efficiency. Check

the correct correspondence of the connectors of the chosen model to those of the machine in use, and proceed to connect the 2 tubes respectively to the pneumatic that drives the movement of the guillotine and to the aspiration of the machine. In some models the pneumatic way will be doubled, while the aspiration way remains unchanged. Connect the connectors carefully and make sure the connection is perfect and complete. Please note that AKTIVE offers a wide range of models with different working pressures and different connections. It is up to the surgeon, and only the surgeon, to decide which one fits the instrument in use, without AKTIVE having any responsibility for this. Depending on the surgical technique the VITRECTOMY CUTTER to reach the area of Vitreous Humor to be removed must be used through an incision made by the surgeon on the eyeball. This incision can be made with a scalpel blade and therefore the VITRECTOMY CUTTER will pass through the scleral tissue to reach the posterior pole in direct contact with the sclera itself. Differently, in the case of minimally invasive surgery using a Trocar, THE VITRECTOMY CUTTER can be introduced directly through the trans-scleral trocar.

DISPOSAL

The used device must be disposed of in accordance with standard hospital procedures to control the risk of contamination and with care taken to protect operators from inadvertent injury during disposal. Secondary packaging, IFUs, blisters and Tyvek® pouches, are medical waste that can be assimilated as municipal waste and therefore disposed of according to procedures applied by the hospital and in accordance with municipal sorting rules.

PACKAGING AND STERILISATION

Single, sterile pack, packaged in aseptically openable blister packs.
ETO Sterilisation

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
Rev. 9 de 20.10.2022

RANGE OF MODELS and WORK SPECIFICS IDENTIFICATION					
PSI	Max CPM	20G	23G	25G	27G
17-18	2.000/4.000 Single/dual cut	nd	1401-2323 Single cut	nd	nd
20	800/1.600 Single/dual cut	AK-0201/AA single cut AK-0201/AA2 dual cut	AK-0231/AA single cut AK-0231/AA2 dual cut	nd	nd
20	800/1.600 Single/dual cut	AK-0215/AB single cut AK-0215/AB2 dual cut	AK-0245/AB single cut AK-0245/AB2 dual cut	AK-0275/AB single cut AK-0275/AB2 dual cut	AK-0294/AB single cut AK-0294/AB2 dual cut
20	2.000/4.000 Single/dual cut	AK-0200/BA single cut AK-0200/BA2 dual cut	AK-0230/BA single cut AK-0230/BA2 dual cut	AK-0260/BA single cut AK-0260/BA2 dual cut	AK-0286/BA single cut AK-0286/BA2 dual cut
25	3.000/6.000 Single/dual cut	AK-0206/AA single cut AK-0206/AA2 dual cut	AK-0236/AA single cut AK-0236/AA2 dual cut	AK-0266/AA single cut AK-0266/AA2 dual cut	AK-0288/AA single cut AK-0288/AA2 dual cut
25	5.500/11.000 Single/dual cut	AK-0208/AB single cut AK-0208/AB2 dual cut	AK-0238/AB single cut AK-0238/AB2 dual cut	AK-0268/AB single cut AK-0268/AB2 dual cut	AK-290/AB single cut AK- 290/AB2 dual cut
26-27	6.000/12.000 Single/dual cut	AK-0207/BA single cut AK-0207/BA2 dual cut	AK-0237/BA single cut AK-0237/BA2 dual cut	AK-0267/BA single cut AK-0267/BA2 dual cut	AK-0289/BA single cut AK-0289/BA2 dual cut
29	2.500/5.000 Single/dual cut	AK-0203/AB single cut AK-0203/AB2 dual cut	AK-0233/AB single cut AK-0233/AB2 dual cut	nd	nd
30	2.000/4.000 Single/dual cut	AK-0210/AB single cut AK-0210/AB2 dual cut	AK-0240/AB single cut AK-0240/AB2 dual cut	AK-0270/AB single cut AK-0270/AB2 dual cut	AK-0292/AB single cut AK-0292/AB2 dual cut
30	5.000/10.000 Single/dual cut	AK-0204/AB single cut AK-0204/AB2 dual cut	AK-0234/AB single cut AK-0234/AB2 dual cut	AK-0264/AB single cut AK-0264/AB2 dual cut	AK-0284/AB single cut AK-0284/AB2 dual cut
33	800/1.600 Single/dual cut	AK-0202/AA single cut AK-0202/AA2 dual cut	AK-0232/AA single cut AK-0232/AA2 dual cut	nd	nd
35-37	5.000/10.000 Single/dual cut	AK-0209/EAA single cut AK-0209/EAA2 dual cut	AK-0239/EAA single cut AK-0239/EAA2 dual cut	AK-0269/EAA single cut AK-0269/EAA2 dual cut	AK-0291/EAA single cut AK-0291/EAA2 dual cut
42	2.000/4.000 Single/dual cut	AK-0220/AC single cut AK-0220/AC2 dual cut	AK-0250/AC single cut AK-0250/AC2 dual cut	AK-0280/AC single cut AK-0280/AC2 dual cut	AK-0295/AC single cut AK-0295/AC2 dual cut
42	5.000/10.000 Single/dual cut	AK-0205/BF single cut AK-0205/BF2 dual cut	AK-0235/BF single cut AK-0235/BF2 dual cut	AK-0265/BF single cut AK- 0265/BF2 dual cut	AK-0287/BF single cut AK-0287/BF2 dual cut
42	6.000/12.000 Single/dual cut	AK-0219/AC single cut AK-0219/AC2 dual cut	AK-0249/AC single cut AK-0249/AC2 dual cut	AK-0279/AC single cut AK-0279/AC2 dual cut	AK-0297/AC single cut AK-0297/AC2 dual cut
42	8.000/16.000 Single/dual cut	AK-0221/GC single cut AK-0221/GC2 dual cut	AK-0251/GC single cut AK-0251/GC2 dual cut	AK-0281/GC single cut AK-0281/GC2 dual cut	AK-0298/GC single cut AK-0298/GC2 dual cut
43	6.000/12.000 Single/dual cut	AK-0211/AB single cut AK-0211/AB2 dual cut	AK-0241/AB single cut AK-0241/AB2 dual cut	AK-0271/AB single cut AK-0271/AB2 dual cut	AK-0293/AB single cut AK-0293/AB2 dual cut
52	800/1.600 Single/dual cut	AK-0226/AD single cut AK-0226/AD2 dual cut	AK-0256/AD single cut AK-0256/AD2 dual cut	nd	nd
52	7.500/15.000 Single/dual cut	AK-0225/ABB single cut AK-0225/ABB2 dual cut	AK-0255/ABB single cut AK-0255/ABB2 dual cut	AK-0285/ABB single cut AK-0285/ABB2 dual cut	AK-0299/ABB single cut AK-0299/ABB2 dual cut

Aktive Srl
V. G. Delitala 106
00173 Roma Italia
www.aktive-corp.com

ZKTVO85