

## MULTI FUNCTIONAL HAND-PIECES

AQUAFEEL	Skin Exfoliation & Moisture
ELECTROPORATION	Transdermal Delivery
MICRO-CURRENT	Facial Lifting
COOLING	Skin Relaxation

## SYSTEM SPECIFICATIONS

Input Voltage	100-240V~, 50/60Hz
Applicator	Aquapeel / Electroporation / Y Lifting / Cryo
Panel	Touch panel
Dimension	364 x 552 x 456mm (WxDxH)
Weight	14kg



AQUAPURE functions include:

- Skin Rejuvenation
- Exfoliation & Extraction
- Deep Pore Cleansing
- Skin Elasticity

AQB022017

- Non-Ablative & Non-Invasive
- Direct Serum Infusion
- LED Color Therapy
- Customized Procedures
- Patented & Proven Formula



Uniform  
Spray  
Control  
Technology



**Cluederm™** LIFTING PATCH



Distribution Switzerland:



www.lasermed.ch

Entire Beauty Solution  
**Cluederm™**

www.cluederm.com

Images and texts are intellectual property of Cluederm. Copying of this material can be subject to charges of both civil and criminal law of legal justice.  
Copyright to Cluederm © 2017

Go Pure · Be Youth

**AQUAPURE**

Smart Genuine Facial Care System

**Cluederm™**  
Entire Beauty Solution

Go Pure · Be Youth

# AQUAPURE

Smart Genuine Facial Care System



## Optimized User-Focused Technology

To address a multitude of skin concerns


### Uniform Spray Control Technology

The AQUAPEEL spray control knob delivers equal distribution of applied solutions at a range of 10 levels. Addressing this concern within a developing microdermabrasion technology maximizes treatment efficiency and results.



Other Models

VS



AQUAPURE



### User-Focused Smart Auto Reaction

The Smart Auto-Reaction Interface detects and enables automatic recognition to all activated hand-pieces.

### Subtle & Stable Electroporation Process

The precision of the ELECTROPORATION hand-piece transfers controlled stimulation to skin without signs of irritation.

### Hand-pieces

· AQUAPEEL



· ELECTROPORATION



· MICRO-CURRENT



· COOLING / HEATING



### Clinical Cases



Baseline

Immediately Post Procedure



Baseline

Immediately Post Procedure



Baseline

Immediately Post Procedure



Baseline

Immediately Post Procedure

### Proven Cosmetic Solution

A combination of natural extracts, Hyaluronic Acid, and a complex of Peptides complement the patented technology, to provide the highest grade of laboriously selected and nutritious ingredients to consumers.



Seaweed Complex



Propolis Extract



Centella Asiatica TECA



Peptide Complex

### AQUASOLUTION Range

PEEL+




The first step of treatment is the process of non-irritating skin scaling which facilitates accelerated exfoliation, the removal of dead skin cells and other impurities.

SEBO+



The second step of treatment is for sebum dissolution, treatment of blackheads, and effective pore control. This assists the skin to balance between oily and moisture textures.

REJUVE+



The final step of treatment promotes skin elasticity by inducing collagen and elastin regeneration and improving the inhibition of moisture loss by improving water retention.

CLEANSYS



To provide a convenient cleaning agent specifically made for the needs of the system, CLEANSYS assists in easy maintenance procedures for the AQUAPURE.