

D-ACTOR® 100 »ultra«



D-ACTOR® mode:

- Radial acoustic wave:
pressure: 0.3 – 5.0 bar_{eff}
frequency: 1 – 21 Hz
(depending on pressure)
- Radial handpiece FALCON® with integrated display and »Skin touch« function for rapid treatment navigation
- Vibration therapy (optional), V-ACTOR®: 1 – 31 Hz
- 10" touch screen (optional)

D-ACTOR® 200 »ultra«



D-ACTOR® mode:

- Radial acoustic wave:
pressure: 0.3 – 5.0 bar_{eff}
frequency: 1 – 21 Hz
- Radial handpiece FALCON® with integrated display and »Skin touch« function for rapid treatment navigation
- Two D-ACTOR® handpiece connectors
- Vibration therapy (optional), V-ACTOR®: 1 – 35 Hz
- Vacuum therapy (optional), VACU-ACTOR® technology
- 10" touch screen (optional)

CELLACTOR® SC1 T-TOP »ultra«



C-ACTOR® mode:

- Focused acoustic wave:
0.03 – 1.24 mJ/mm²
- Focused handpiece C-ACTOR® SEPIA® with integrated display for rapid treatment navigation
- 10" touch screen (optional)

CELLACTOR® SC1 TOWER »ultra«



C-ACTOR® mode:

- See CELLACTOR® SC1 T-TOP »ultra«
- D-ACTOR® mode:**
- See D-ACTOR® 200 »ultra«
 - Integrated ultrasound imaging, Colour Doppler (optional)
 - Drawer module »Store Case« (optional)
 - 15,4" touch screen