



ADVANCED  
TECHNOLOGY FOR  
FRACTIONAL SKIN  
REJUVENATION



## WHAT IS MEZOTIX™?

- A novel non ablative fractional skin rejuvenation system powered by Novoxel's patented titanium tip technology
- Mezotix treats delicate facial skin, including peri-orbital, eyelids, the neck and décolleté safely and quickly with low pain
- Employs pure heat and is radiation-free

## OUR BENEFITS



### MULTIPLE TREATMENTS

Skin rejuvenation, crow's feet, open pores, acne, stretch marks ...



### MULTIPLE AREAS

Face, neck, décolleté, hands, arms & body



### LOW PAIN

Analgesic creams are not required



### STERILE TITANIUM TIP

The tip is biocompatible and sterile



### SMALL & PORTABLE

Small, light, easy to carry



### NO RADIATION

Mezotix employs pure natural heat

> · Glabellar  
"frown" lines

> · Eyelid lines

> · Acne scars

> · Creased chin

> · Décolleté

> · Forehead  
"worry" lines

> · Crow's feet  
(Peri-orbital lines)

> · Bunny lines

> · Nasal labial folds  
(Smile lines)

> · Lip lines

> · Neck folds





# TREATMENTS

- Spacing treatment sessions - 1 month apart
- Treatment Program: 3 - 4 sessions
- After treatment skin is slightly swollen and red for 3 to 4 hours
- Patients experience erythema and heat for 2 to 3 hours
- There is typically no downtime
- Crusting appears the next day, disappears after about 3 days
- Improvement is usually visible after 1<sup>st</sup> month, 1<sup>st</sup> session
  - Improvement in skin texture, younger look, wrinkle reduction
  - Better definition of eyelids
  - Patients express satisfaction









**CORINA MÜLLER,  
JE'SCAVA**

(NUREMBERG, GERMANY)

## CUSTOMER FEEDBACK

"Mezotix is a very compact and versatile device. It's quiet and easy to use. Treatment with Mezotix is extremely fast, very safe and results are outstanding. For me Mezotix is the jewel of anti-aging procedures. My patients appreciate their results and have greater comfort, feeling no pain. They resume their daily social activity quickly. The combination of economic advantages and versatility of the Mezotix enables a quick ROI."

## TESTIMONIALS

“ The Mezotix treatment was without any pain or discomfort and I was able to work on the following day. ”

“ I just experienced a warm sensation in the treatment area similar to a sunburn which dissipates after about an hour. ”

“ My Mezotix treatment was very quick. It takes just about 15 minutes for my hands, face and décolleté. ”

“ I tried so many different treatments but with Mezotix the result was outstanding like no other before. ”

## BETTER FOR CLINIC


- Excellent clinical results
- Enhanced client experience
- Faster treatment time
- Easy to operate
- Treat darker skin tones

## BETTER FOR PATIENT

- Less pain
- Less charring & crusting
- Less downtime & post-care
- Faster treatment time
- Safer

# TECHNICAL SPECIFICATIONS

Energy Type:	Direct heating
Tip Material:	Titanium
Tip Coverage Area:	1 cm <sup>2</sup> Standard Tip, 0.3 cm <sup>2</sup> Periorbital
System Weight (w/o HP):	4.5 kg
Handpiece Weight:	330 grams
System Dimensions (w/o HP):	L 26 x W 24 x H 24 (cm)
Operating Voltage:	100-240 VAC ±10% 5A 50/60Hz



## TITANIUM TIP PROPERTIES

- >• DURABLE
- >• STERILE
- >• CLEAN
- >• PRECISE



Gatower Str. 124

13595 Berlin, Germany

Tel: +49 (0) 30.548.11195

Fax: +49 (0) 30.548.16951

Email: [Info@novoxel.com](mailto:Info@novoxel.com)

## Customer Support

Tel: +49 (0) 30.548.11195

Email: [supportGmbH@novoxel.com](mailto:supportGmbH@novoxel.com)

Novoxel GmbH. All rights reserved. Novoxel GmbH., its logo Novoxel®, is a trademarks of Novoxel.  
Product specifications are subject to change without notice.

\*Patent pending.

### Distribution Switzerland:



LASERMED AG  
Roggwil TG  
Corminboeuf FR  
[info@lasermed.ch](mailto:info@lasermed.ch)  
[www.lasermed.ch](http://www.lasermed.ch)

P.N. 1284 / Brochure Mezotix / Version 01