

Esophageal Manometry System



ManoScan[™] High Resolution Manometry Diagnosing with definition

ManoScan[™] ESO

ManoScan ESO provides a complete physiological mapping of esophageal motor function, from the pharynx to the stomach, with a single placement of a catheter. This advanced diagnostic technology allows physicians to better diagnose disease states such as dysphagia, achalasia, and hiatal hernia. The procedure is easier for the clinician to perform and is more patient-friendly than conventional manometry.

- The only system with published normative clinical data²
- HRM can precisely quantify the contractions of the esophagus and its sphincters⁴
- Most studies completed in 10 minutes or less and require minimal specialized training³

 2 Pandolfino JE et al Neurogastroenterol Motil 2009;21:796-806.

³ Kahrilas PJ et al. Am J Gastroenterol 2010;105:981-987.

⁴ Bansal A et al. Curr Opin Gastroenterol 2010;26;344-351.













Achalasia Type II

ManoScan[™] ESO Z

ManoScan ESO Z provides circumferential assessment of bolus movement as well as physiological mapping of the esophageal motor function, from the pharynx to the stomach, with a single placement of the catheter.

- The incorporation of impedance measurement with HRM maps improves the ability to predict the success or failure of bolus movement through the esophagus
- This technology aids physicians in better understanding the causes of dysmotility, such as achalasia, dysphagia, and reflux



Normal Swallow with Impedance

ManoScan[™] V

The ManoScan video module works in conjunction with high resolution manometry to allow for synchronized, simultaneous video and pressure collection providing a previously unseen diagnostic picture. When used with ManoScan ESO, this module pairs pressure mapping with real-time video visualization of swallow coordination.

- Fluoroscopic studies can provide complementary information to HRM in order to confirm diagnosis and treatment
- Provides tremendous potential for pharyngeal biofeedback retraining in stroke victims and cancer patients



Full Featured Workstation

- Portable trolley system
- LCD flat panel touchscreen with articulating arm
- Modular data acquisition controller
- Windows[®]-based operating system
- LAN connection and WiFi enabled
- Integrated catheter auto-calibration system
- Built-in storage drawer
- Large lockable wheels
- Patient isolation transformer
- High speed quality printer

ManoView[™] Software

The ManoView[™] software provides an intuitive suite of manometry study tools enabling physicians to effectively identify motility disorders.

- Advanced tools yield precise measurement and comprehensive data analysis
- Anatomical profile display includes graphical pointers to identify landmarks including LES, UES, and PIP
- eSleeve[™] function instantly measures and ensures sphincter barrier pressures are correctly recorded despite movement of the LES/EGJ during swallowing
- High resolution and conventional displays provide versatile and complete motility visualization
- ManoView software can be installed on any Windows[®]based computer enabling clinicians to review studies remotely







ManoScan Catheter



ManoScan[™] HRM Catheters

ManoScan HRM catheters incorporate the very latest advancements in sensing technology.

- ManoScan HRM catheters work with ManoShield disposable sanitary sheath to minimize the risk of cross-contamination
- With 36 channels providing 432 points of measurement, the ManoScan ESO catheter provides the highest resolution of any available manometry catheter

ManoScan[™] Z Catheters

ManoScan Z catheters incorporate the latest advancements in sensing technology.

- All sensors are true circumferential
- 36 pressure channels spaced 1 cm apart, create a pressure image from pharynx to stomach
- 18 impedance channels display bolus transition from pharynx to stomach

ManoScan Z Catheter

ManoShield[™] Disposable Catheter Sheath

Micro-thin disposable sheath design creates a protective barrier between the catheter or probe and the patient to enhance patient safety during HRM procedures.

- Creates more efficient workflow and minimizes catheter deterioration by significantly reducing need for cleaning and disinfecting between procedures
- Minimizes risk of cross-contamination, helping clinicians comply with health care patient safety requirements
- All materials are latex-free and hypo-allergenic



ManoShield

Diagnosing with Definition



Jeffrey H. Peters, MD Chairman, Department of Surgery University of Rochester, New York



Copyright ©2001-2011 Given Imaging Ltd. GIVEN, GIVEN & Design, PILLCAM, PILLCAM & Logo, PILLCAM IMAGING CAPSULE & Design, PILLCAM EXPRESS, AGILE, RAPID, RAPID ACCESS, ORDERWIN, FINGERS HOLDING A CAPSULE & Logo, FINGERS HOLDING PILLCAM CAPSULE & Logo, ICCE, ICCE Logos, International Conference on Capsule Endoscopy, BRAVO, BRAVO PH SYSTEM, ENDONETICS, VUESPAN, VERSAFLEX, GEROFLEX, REPHLUX TRACER, GASTROTRAC, BILITEC, DIGITRAPPER, PHERSAFLEX, MANOSCAN, MANOSCAN 360, MANOSCAN Z, MANOSCAN 3D, MANOSCAN HD, MANOSCAN V, MANOSHIELD, MANOSHIELD AR, MANOVIEW, MANOVIEW Z, ACCUTRAC, ACCUTRAC Z, ACCUVIEW, ACCUVIEW Z, ION, INSERTASSIST, BOLUSVIEW, and POLYGRAF ID are Trademarks and/or Registered Trademarks of Given Imaging Ltd., its subsidiaries and/or affiliates in the United States and/or other countries. All other company or product names are the trademarks or registered trademarks of their respective holders. All rights not expressly granted are reserved.

Shares of Given Imaging trade on the NASDAQ Global Market and on the Tel-Aviv Stock Exchange under the symbol "GIVN".



www.givenimaging.com